




# Compact Q

Real-Time PCR Detection System




Powered by Compact Q,  
Process 16 samples  
within 40 minutes while  
maintaining the gold  
standard accuracy.


# Product Performance

**Experiment 1:** Experiment of distinguishing  
between 2000 copy concentration and  
1000 copy concentration

**Experiment 2:** Standard curve  
performance test of HBV of 5μl



Experimental results show that the Ct difference is extremely close to the perfect theoretical value 1, Compact Q has accurate detection performance, and can perfectly distinguish the detection data at different concentrations.



Experimental results show that Compact Q has excellent concentration gradient detection performance and good curve fitting effect.

# Key Features & Benefits



Advanced scanning technology  
with excellent SNR



Exquisite and compact exterior  
with strong and durable structure



Cutting-edge detection technology  
and most stable signal



Touch operation interface,  
easy to operate




Automatic gain function,  
accurate and reliable data



Complete function configuration  
& stable equipment performance

Write comes with  
MylabCares auto  
reporting software

# Machine Specifications



Product Name	Compact Q-Fluorescence Quantitative PCR Detection System			
Product model	FQD-16B			
Sample capacity	16x0.2ml single tube (transparent tube)		0.2ml 8 strip tube (transparent tube)	
Reaction system	5-100μl			
Dynamics range	1-10 <sup>10</sup> Copies/L			
Excitation wavelength	400~700nm			
Fluorescence Emission wavelength	450~750nm			
Temperature range	0-100° C (Minimum increment: 0.1°C			
Temperature accuracy	± 0.1°C			
Temperature display resolution	0.1°C			
Temperature control accuracy	± 0.1°C			
Temperature uniformity	≤ ± 0.25°C (95°C)			
Maximum temperature rate	5.0°C/s			
Temperature control mode	Block model 、 Tube model			
Hot-lid temperature range	30-110°C (adjustable)			
Detection path	F1	F2	F3	F4
Fluorescence dyes	FAM 、 SYBR GreenI...	VIC 、 HEX 、 TET 、 JOE 、 TAMRA 、 CY3 、 NED...	ROX 、 TEXAS-RED...	CY5 、 Quasar-670...
Special functions	Qualitative/absolute quantitative, relative quantitative, SP analysis, fusion curve method genotyping, standard curve, fusion curve, allele identification, gradient, HRM, automatic gain, user defined parameters, etc			
Input power	100-240V, 50/60Hz, (A power adapter is required 180w/24V)			
Signal interface	Netport (connected with PC)			
Dimension(mm)	240×180×160 (L×W×H)			
Net weight	4Kg			
Safeguard and alarm	Support temperature sensor short circuit, open circuit & over-temperature protection & alarm			
Certification	CE, In line with the ROHS requirements			

**Mylab Discovery Solutions Pvt. Ltd.**

Plot No. 99-B, Lonavala, Industrial Co-operative Estate Ltd., Nangargaon, Lonavala, Pune, Maharashtra 410401, INDIA.

**Email:** [info@mylabglobal.com](mailto:info@mylabglobal.com) | **Web:** [www.mylabglobal.com](http://www.mylabglobal.com)

**For more information  
on the products contact:**

**Varad Gawde**  
**+91 98507 82682**  
[Varad@mylabglobal.com](mailto:Varad@mylabglobal.com)

**Mohit Agrawal**  
**+91 99303 75502**  
[Mohit@mylabglobal.com](mailto:Mohit@mylabglobal.com)

**Karthik Balasubramanian**  
**+91 91675 90922**  
[Karthik.b@mylabglobal.com](mailto:Karthik.b@mylabglobal.com)

**CoviSelf** | **CoviSwift** | **PathoDetect** | **NATSpert** | **PathoCatch** | **OncoScreen** | **DiscoverSeries** | **AgroDex** | **VetScreen** | **CompactXL**

[mylabglobal.com](http://mylabglobal.com)



[mylabsolutions](https://twitter.com/mylabsolutions)



[mylabsolutions](https://www.instagram.com/mylabsolutions)



[mylabdiscoverysolutions](https://www.facebook.com/mylabdiscoverysolutions)



[mylabdiscoverysolutions](https://www.linkedin.com/company/mylabdiscoverysolutions)

**Disclaimer:** All rights reserved. All content (texts, trademarks, illustrations, photos, graphics, files, designs, arrangements, instruments etc.) are protected by our Intellectual Property (IP) including Patents, Copyright, Trademark and other protective laws. Any attempt for infringement of IP vested in our products, for any purpose whatsoever, without taking prior written license from Mylab would amount to infringement and may be harmful involving high risk to life and may result in fatal injury. **For enquiries related to IP licensing, please write to [legal@mylabglobal.com](mailto:legal@mylabglobal.com)**